

THE SHELVES CAN ALSO BE ORDERED WITH AN ANTIMICROBIAL SUPPLEMENT!

Upgrade around your sterile storage with professional antimicrobial storage units!



WEIGHT CAPACITY UP TO 220LBS (100KG) PER SHELF



removes + 99% of bacteria in 24 hours



Continuous protection for the entire life of the shelving system



Temperature range: from -22°F (-30°C) to +194°F (+90°C)



Impenetrable defense against microbes, bacteria, fungi and molds



treatment suitable for direct food contact




WILL NOT RUST
Anodized Aluminum Structure

SHELVES COME IN THREE DIMENSIONS

SHELVES

DESCRIPTION

ITEM NUMBER




Hygieia Shelf UU

HYG-SHELF-AM-UU

1 Pcs/bx

HDPE (RAL 6034) & Anodized Aluminum Rails

* Requires one LS set of 3 per rack section: HYG-LSU-UU




Hygieia Shelf UUU

HYG-SHELF-AM-UUU

1 Pcs/bx

HDPE (RAL 6034) & Anodized Aluminum Rails

* Requires one LS set of 3 per rack section: HYG-LSU-UUU



Hygieia Shelf EE

HYG-SHELF-AM-EE

1 Pcs/bx

HDPE (RAL 6034) & Anodized Aluminum Rails

* Requires one LS set of 3 per rack section: HYG-LSU-EE

ANTIMICROBIAL MODULES

INFECTION CONTROL storage (ABS with antibacterial agents):




Tray 10cm (4\") Deep

T104060ND-AM



Tray 10cm (4\") Deep

T104060-AM



Tray 20cm (8\") Deep

T204060-AM