

INTEGRATED HEAVY-DUTY PLASTIC SHELVING FOR THE HYGIEIA ALUMINUM RACK SYSTEM

Designed for sterile locations including the OR and SPD standard and antimicrobial shelves available



WEIGHT CAPACITY UP TO 220LBS (100KG) PER SHELF



Lightweight



Durable



No special tools
needed for assembly



Easy to clean



Temperature range:
from -22°F (-30°C) to
+194°F (+90°C)




WILL NOT RUST
Anodized Aluminum
Structure

SHELVES COME IN THREE DIMENSIONS

SHELVES

DESCRIPTION

ITEM NUMBER




Hygieia Shelf UU

HYG-SHELF-UU

1 Pcs/bx

HDPE (RAL 7040) &
Anodized Aluminum Rails

* Requires one LS set
of 3 per rack section:
HYG-LSU-UU




Hygieia Shelf UUU

HYG-SHELF-UUU

1 Pcs/bx

HDPE (RAL 7040) &
Anodized Aluminum Rails

* Requires one LS set
of 3 per rack section:
HYG-LSU-UUU



Hygieia Shelf EE

HYG-SHELF-EE

1 Pcs/bx

HDPE (RAL 7040) &
Anodized Aluminum Rails

* Requires one LS set
of 3 per rack section:
HYG-LSU-EE

LATERAL SUPPORT


Hygieia lateRAL support UU Anodized aluminum



HYG-LSU-UU

1 set (3 rails/bx)

Hygieia lateRAL support UUU Anodized aluminum



HYG-LSU-UUU

1 set (3 rails/bx)

Hygieia lateRAL support EE Anodized aluminum



HYG-LSU-EE

1 set (3 rails/bx)