EMERGENCY/ CRASH CARTS COMPLETE SOLUTION



ELEVATING EMERGENCY RESPONSE
WITH ADAPTIVE, HYGIENE-CENTRIC CRASH CARTS



Stay updated about your Crash Carts with Audit notifications.

See pg 11 for more info



OPTIMAL INNOVATION

With over 20 years dedicated to studying and understanding the needs of healthcare professionals, our solutions are born from a wealth of expertise. The user's needs serve as our guiding light, shaping innovations that empower and deliver excellence. From swift and lightweight to resilient and premium quality, we offer everything essential for emergencies. Join us at Pegasus Medical in our journey to create a more efficient and fast emergency response!

WORKFLOW



Contact us

info@pegasusmedical.net 1-888-276-4750



Configurate

You can choose accessories, sizes, colors



Approval

Review our offer & give your approval



Ready to ship

We handle installation and delivery processes



Your work optimized, efficient, and safe



READY TO SHIP FOR ALL THE WORLD

REVOLUTIONARY CRASH CARTS

Diverse range of cutting-edge products, each designed to meet specific needs.



POLY LINE

100% antimicrobial polymer

ECOLINE 40by40

Perfect for narrow spaces!

UNIFIED EMERGENCY SOLUTIONS

We cater to all emergency scenarios with our diverse range of specialized emergency carts, tailored to every situation and need

POLY PEDIATRIC

MULTIPLE CARE & SECURE

Individual locks allow for simultaneous use in multiple emergencies and enable easy and fast restocking



PEDIATRIC

ADULT



CCANO

SCAN & AUDIT

Pg 11 for more info





ECOLINE 40by40

✓ ADAPTABLE TO LIMITED SPACES

Small footprint allows it to navigate tight spaces



POLY EMERGENCY

J QUICK OPENING

With one pull you can open all the drawers at once

✓ FAST & LIGHTWEIGHT

Lightweight and easy to maneuver, made from antimicrobial polymers



POLY CART

Innovating for the future in emergency care, our line now introduces a groundbreaking 100% antimicrobial solution — a revolution in all medical facilities

✓ SECURE & LABELED

The cart features two different types of locks — an external lock and a seal lock — with breakaway seals and a drawer labeling system, ensuring fast identification and usage of the cart's content

LIGHTWEIGHT & RESISTENCE

Crafted from advanced polymer materials, this cart not only prioritizes durability and resistance but also ensures exceptional lightweight performance

OVERIDE LOCK

Disable all locks in case of malfunction ensuring access to the cart







Poly Cart with Seal Lock



Poly Cart with External Lock

DRAWER LABELING SYSTEM



POLY CART

Enhanced organization, facilitates inventory management and reduces response time



CLEAR & ORGANIZED

Equipped internally with baskets of various dimensions and multiple dividers, the cart incorporates an internal labeling system





FAST & SECURE

The emergency hospital cart features four 5-inch swivel wheels for smooth rotation and quick movement, with two front wheels providing stability during stationary periods and critical procedures. Optionally, it can be equipped with a directional lock for fast, straight movement in long corridors



TAILORED SOLUTION

Tailor our emergency cart to your exact specifications with customizable features including height, colors, locks, and a range of accessories. This ensures a uniquely crafted solution that meets the specific needs of your healthcare environment



COLORS







PLY46-444449-RR 46" | 117 cm



PLY41-44449-RR 41" | 104 cm



PLY36-4449-RR 36" 92 cm




84 cm x 58 cm (LxW) 33" x 23" (LxW) THREE HEIGHTS

ACCESSORIES

A range of accessories to fit your exigencies



IV pole PLY-IVPOLE





Wheels



Swivel 5" castor with directional lock (Optional) PLY-S5CDT



Swivel 4" dual castor total lock & directional lock (Optional - cart height decreases by 1") PLY-D4CST-DT

PEDIATRIC POLY CART

Optimizing pediatric emergency response



Individual locks allow for simultaneous use in multiple emergencies and enable easy and fast restocking


✓ QUICK & CLEAR

Fast opening with one pull, and easy to see which drawers were utilized


✓ RELIABLE

Being independent of an internal mechanical locking system guarantees lifetime access





CRASH CART CONTROL



DECREASE LIABILITY

Every audit is logged with team member details and can be filtered by relevant criteria, effectively reducing hospital liability

/ AUDIT **NOTIFICATIONS**

Receive alerts regarding failure to perform daily/monthly cart audits, eliminating human errors and forgetfulness



EXPIRATION PUSH ALERTS

Receive notifications for expiring or low supplies, and stay informed about the status and functionality of all equipment



demo@intraposition.com Password: Demo1234!





UNIQUE ID

Every crash cart has its unique QR code that provides online access to all reports and audits for that specific cart



SIMPLICITY

Scan the code for instant access to relevant audit forms, optimizing staff time management



GREEN INITIATIVE

Decrease your carbon footprint by going paperless


ECOLINE 40by 40

✓ SMALL FOOTPRINT, **BIG IMPACT**


Compact design allows it to navigate tight spaces effortlessly, making it the perfect fit for any healthcare environment

/ INTEGRATED POLES


We've seamlessly integrated poles for top accessories into the frame, enhancing its strentgh and stability







ACCESSORIES







PEGASUS MEDICAL CONCEPTS

905 E Rose St, Lakeland, Florida 33801, US

TEL. 888-276-4750

info@pegasusmedical.net

WWW.PEGASUSMEDICAL.NET