

NO *Touch*™



THIS


NOT THIS

**CONTAMINATION AND TEARS
ARE A THING OF THE PAST**

NO Touch™

CONTAMINATION AND TEARS ARE A THING OF THE PAST WITH NO-TOUCH™... FROM THE SPD TO THE OR, WITH VIRTUALLY NO HUMAN INTERACTION

The No Touch system eliminates human interaction during instrument reprocessing thanks to our dedicated ISO Standard Stainless Steel shelves and baskets and our functional lines of open and closed transportations carts.



STEP 1
Clean Instrument Wrapping



STEP 2
Wrapped Instruments are placed on No Touch Shelf and into Sterilizer (shelf is handled, pack is not touched)



STEP 3
No Touch Shelf (w/sterile pack) is transported to Sterile Storage (shelf is handled, pack is not touched)



STEP 4
No Touch Shelf (w/sterile pack) is transported to OR to be Opened by Nurse


Mobile & Stationary Storage Systems



Stainless steel Transport Cart



Open Transportation Trolley



Sterile Instrument Rack Storage

STAINLESS STEEL MODULES


Modules suitable for sterilization




60X40cm (24"X16")
10cm (4") Deep



60X40cm (24"X16")
20cm (8") Deep



Flat 60X40cm (24"X16")



Flat 30X40cm (12"X16")



THIS

Eliminate contamination and tears during transportation by touching only the carrier shelf, and not the sterile pack



NOT THIS

Maximum sterility during transportation

Safeguards patients from infection

NO *Touch*™



CONTACT US FOR MORE DETAILS

+1 888-276-4750 | info@pegasusmedical.net | www.pegasus-notouch.com

USA / CANADA
905 E. Rose St Lakeland
FL 33801 - USA

EUROPE / MIDDLE EAST / ASIA
Damzigt 13, 3454 PS De Meern
The Netherlands