

Fit your Ecoline crash cart with broselow color coding

The set includes 8 handles each in a different color and complies with the Broselow standard

The new Pegasus Handle Set is designed to give you the ability to fit a basic Ecoline cart into a pediatric emergency cart by simply replacing the handles. The cart recommended configuration should include 9 drawers*

Item number: **CRTA-MKHND-BRS-SET**



THE HANDLE SET INCLUDES THE FOLLOWING COLORS

- | | |
|------------|---------------|
| 1 X Pink | 1 X White |
| 1 X Red | 1 X Dark blue |
| 1 X Purple | 1 X Orange |
| 1 X Yellow | 1 X Green |

*Use one light gray handle that comes by default in each Ecoline cart for the top drawer to complete the Broselow standard.