


APOLLO MEDICATION CART



The Apollo Medication Cart offers an ideal solution for medication distribution at the Point of Care. The cart has an extremely light structure, thanks to its combination of an aluminum skeleton with polymeric walls

COMPUTER ARM

Can carry up to 18kg, with tilt position adjustment



KEYBOARD DRAWER

The drawer can be furnished with a 10cm tray or ergonomic keyboard inlay



MEDICATION CASSETTE

The cart is equipped with light-weight ergonomic removable medication cassette, mounted on a set of full-length retracting telescopic slides



The cassette can be equipped with 12 to 18 dividable medication bins. Each bin can be labeled and equipped with a label cover, so no need for sticky labels

STORAGE DRAWER

The lower drawer is equipped with a 20cm tray for general storage or liquid/ bottle storage

The Pegasus modular medication cassette has a footprint of 60X40cm therefore it can easily slide into a transportation cart that are designed to accommodate ISO 3394 modular systems



FOR MORE DETAILS

+1 888-276-4750 | info@pegasusmedical.net | www.pegasusmedical.net