




- ALL PLASTIC BOX, ABS INJECTED
- EASY TO CLEAN AND WASH
- CLEAR MARKING OF GLOVE SIZE
- MULTI-MOUNTING OPTIONS
- FITS LARGE VARIETY OF BOX SIZES
- STABLE FIXATION OF THE BOX, REAR SPRINGS PUSHING THE BOX INTO POSITION
- QUICK & COMFORTABLE BOX REPLACEMENT



S-GBH-ABS

INTEGRAL SPRINGS


Integral springs pushing the box towards the front of the holder.
Suitable for different glove box sizes.
Maximum box size: 14.5X24.5X8cm



DIN RAIL HOOKS

DIN-HKS-GBH-SET

A set of 2 DIN rail hanging hooks for glove box holder



DIRECT WALL MOUNT

S-GBH-ABS

Horizontally or vertically mounting (screw anchors are not supplied)



PATIENT BED

BED-HKS-GBH-SET

A set of 2 bed hooks for glove box holder



TR-GBH-ABS

Triple box set