


Protect your surgical tools during transportation to/from the procedure rooms to the central sterilization unit

Use the red seal to mark contaminated tools during transportation

L4060-HR


T104060-ND-HR



AVAILABLE OPTIONS


TRAY - # T104060-ND-HR + LID - # L4060-HR

High resistance
40*60cm tray, 10cm high
with lid



TRAY - # T103040-HR + LID - # L3040-HR

High resistance
30*40cm tray, 10cm high
with lid



PSSEAL-R
Disposable Seal, 100 PCS per pack



CONTACT US FOR MORE DETAILS

STAINLESS STEEL MODULES



A dedicated line of high quality stainless-steel baskets and shelves, used in a variety of applications around the central sterilization department and OR suites

Baskets can be nestable/stackable, by rotating the baskets in 180 degrees:




AVAILABLE OPTIONS

ISO Dimensions




60X40cm (24"X16") | 10cm (4") Deep
WS-B104060-SS-WH



Flat 60X40cm (24"X16")
SSS-6040WH




60X40cm (24"X16") | 20cm (8") Deep
WS-B204060-SS




Flat 30X40cm (12"X16")
SSS-3040WH

Related Accessories




Silicone Mesh Mat 42X25cm (16.5"X9.8")
SIMAT-250X420




Silicone Mesh Mat 25X25.5cm (9.8"X10")
SIMAT-250X255

None ISO Dimensions



60X59cm (23.6"X23")
SSS-6059WH



60X29.5cm (23.6"X11.6")
SSS-6029WH

CONTACT US FOR MORE DETAILS

+1 888-276-4750 | info@pegasusmedical.net | www.pegasusmedical.net