

# **NO** *Touch*<sup>TM</sup>



**THIS**


**NOT THIS**

**CONTAMINATION AND TEARS  
OF SURGICAL INSTRUMENT PACKS  
ARE A THING OF THE PAST**

# **NO** *Touch*™

CONTAMINATION AND TEARS ARE A THING OF THE PAST WITH NO-TOUCH™... FROM THE SPD TO THE OR, WITH VIRTUALLY NO HUMAN INTERACTION

The No Touch system eliminates human interaction during instrument transport thanks to our dedicated ISO Standard Stainless Steel shelves and baskets and our functional lines of open and closed transportations carts.



### STEP 1

Clean Instrument Wrapping



### STEP 2

Wrapped Instruments are placed on No Touch Shelf and into Sterilizer (shelf is handled, pack is not touched)



### STEP 3

No Touch Shelf (w/sterile pack) is transported to Sterile Storage (shelf is handled, pack is not touched)



### STEP 4

No Touch Shelf (w/sterile pack) is transported to OR to be Opened by Nurse


# Mobile & Stationary Storage Systems



Aluminum Transport Cart



Open Transportation Trolley



Sterile Instrument Pack Storage

Eliminate contamination and tears during transportation by touching only the carrier shelf, and not the sterile pack

- ✓ Maximum sterility during transportation
- ✓ Safeguards patients from infection



**THIS**



**NOT THIS**

## STAINLESS STEEL MODULES


Modules suitable for sterilization




60X40cm (24"X16") | 10cm (4") Deep  
WS-B104060-SS



60X40cm (24"X16") | 20cm (8") Deep  
WS-B204060-SS



Flat 60X40cm (24"X16")  
SSS-6040WH



Flat 30X40cm (12"X16")  
SSS-3040WH

# NO Touch™



SCAN TO REACH  
THE NO-TOUCH  
PAGE ON OUR SITE



## CONTACT US FOR MORE DETAILS

+1 888-276-4750 | [info@pegasusmedical.net](mailto:info@pegasusmedical.net) | [www.pegasusmedical.net](http://www.pegasusmedical.net)

USA / CANADA  
905 E. Rose St Lakeland  
FL 33801 - USA

EUROPE / MIDDLE EAST / ASIA  
V. Croda Granda 13a 31027  
Spresiano (TV) - Italy